

Gastroenterology

Endoscopy

BIOPSY FORCEPS RE-USABLE

PRODUCT FEATURES

Wide Opening

Large Cup Size

Cups with Back Eye

Grinded Distal Spiral End

Auto - clavable

Smooth Operation Even in Loops

Easy Grip Handle



Cat No.	Deascription	OD	Length
EBF-R118120	Oval cup with teeth	1.8mm	120cm
EBF-R218120	Oval cup with teeth and Needle	1.8mm	120cm
EBF-R123180	Oval cup with teeth	2.3mm	180cm
EBF-R223180	Oval cup with teeth and Needle	2.3mm	180cm
EBF-R123230	Oval cup with teeth	2.3mm	180cm
EBF-R223230	Oval cup with teeth and Needle	2.3mm	180cm

FOREIGN BODY FORCEPS RE-SUSABLE

PRODUCT FEATURES

- Wide Opening
- Grinded Distal Spiral End
- Extra Strong Grip
- Auto-clavable
- Smooth Operation Even in Loops
- Easy Grip Handle



Cat No.	Deascription	OD	Length
EFB - R118120	Rat Tooth Type	1.8mm	120cm
EFB - R123180	Rat Tooth Type	2.3mm	180cm
EFB - R123230	Rat Tooth Type	2.3mm	230cm
EFB - R218120	Alligator Type	1.8 mm	120cm
EFB - R223180	Alligator Type	2.3mm	180cm
EFB - R223230	Alligator Type	2.3mm	230cm
EFB - R618120	Rat Tooth Alligator Type	1.8 mm	120 cm
EFB - R623180	Rat Tooth Alligator Type	2.3mm	180 cm
EFB - R623230	Rat Tooth Alligator Type	2.3mm	230cm

BIOPSY FORCEPS SINGLE USE



PRODUCT FEATURES

Wide Opening

Large Cup Size

Cups with Back Eye

With PTFE Sheath covering

Sterile Single Use

Smooth Operation Even in Loops

Easy Grip Handle

Cat No.	Deascription	OD	Length
EBF - S118120	Toothed Cup	1.8mm	120cm
EBF - S118180	Toothed Cup	1.8mm	180cm
EBF - S123120	Toothed Cup	2.3mm	120cm
EBF - S123180	Toothed Cup	2.3mm	180cm
EBF - S123230	Toothed Cup	2.3mm	230cm
EBF - S223120	Toothed Cup with needle	2.3mm	120cm
EBF - S223180	Toothed Cup with needle	2.3mm	180cm
EBF - S223230	Toothed Cup with needle	2.3mm	230cm

Other sizes are available on request

Rota - Net

FOREIGN BODY NET BASKETS SINGLE USE

PRODUCT FEATURES


Rotatable

Wide Opening

Extra Strong Net

Smooth Operation Even in Loops

Easy Grip Handle



Cat No.	Deascription	Sheath Size	Length
EFB-0730230	Net Type 30 mm	2.3mm	230cm
EFB-0740230	Net Type 40 mm	2.3mm	230cm

POLYPECTOMY SNARES


PRODUCT FEATURES

Braided Wire

Adjustable Size

Easily Rotatable

360 ° Rotatable



Cat No.	Type	Snar size	Sheath size
EPS - 0125230	Oval Shape	25 mm	2.4 mm x 230 cm
EPS - 0130230	Oval Shape	30 mm	2.4 mm x 230 cm
EBS - 0140230	Oval Shape	40 mm	2.4 mm x 230 cm
EBS - 0150230	Oval Shape	50 mm	2.4 mm x 230 cm
EPS - 0230230	Hexagonal Shape	30 mm	2.4 mm x 230 cm

HEMOCLIP SINGLE USE

G-cLip

PRODUCT FEATURES

Pre - loaded clip


360 Rotatable

Wide Opening

Repeated open & close possible

Smooth Operation at difficult location

Easy Grip Handle



Cat No.	Description	Opening	OD	Length
EHS - 0123230	Pre - loaded	11 mm	2.3 mm	230cm
EHS - 0223230	Pre - loaded	13 mm	2.3 mm	230cm
EHS - 0323230	Pre - loaded	16 mm	2.3 mm	230cm

ENDOSCOPY SCLERO THERAPY NEEDLE V - JECTOR

PRODUCT FEATURES

Lock Mechanism

Easy Grip Handle

45 ° Beveled Stainless Steel Needle

Metal Tip

Smooth operation even at difficult

Locations

PTFE inner & outer Tubes



Cat No.	Needle Gauge	Needle Length	Sheath OD	Sheath Length
ESN - 0121180	21G	5 mm	2.3 mm	180 cm
ESN - 0123180	23G	5 mm	2.3 mm	180 cm
ESN - 0121230	21G	5 mm	2.3 mm	230cm
ESN - 0123230	23G	5 mm	2.3 mm	230cm

MULTI BAND LIGATOR

PRODUCT FEATURES

Maximum Endoscopic View


Red Color Trail Band

Penultimate Blue Color Band

Six Usable Bands

Easy Irrigation / Flushing Port

Click Sound With Every Band Release



Cat No.	Scope OD	Band per Barrel	Inclusion
EBL - 0106085-7 (Full set)	9 mm to 11.5 mm	Six + One (Red)	Handle + Barrel + Loading Catheter + Flushing Tube
EBL - 0206085-7 (Mini Pack)	9 mm to 11.5 mm	Six + One (Red)	Barrel + Loading Catheter

Vari - fix[®] Ultra

MULTI BAND LIGATOR

PRODUCT FEATURES

Seven Usable Bands

Extremely Fine Shaped Bands

Different Color Penultimate Band

Easy Irrigation / Flushing Port

Micro Beads Based Thread

Very Thin Natural Yarn



Cat No.	Scope OD	Band per Barrel	Inclusion
EBL - 0107085 (Full Set)	9 mm to 11.5 mm	Seven	Handle + Barrel + Loading Catheter + Flushing Tube
EBL - 0207085 (Mini Set)	9 mm to 11.5 mm	Seven	Barrel + Loading Catheter

DILATION BALLOON

Esophageal , Pyloric / Colonic , Biliary , Achalasia



PRODUCT FEATURES

- 3 stage Dilatation in one Balloon
- Transparent Material
- Metal Marks For Easy Positioning
- Guide Wire Compatible
- Non - Compliant Balloon
- Fast Inflation & Deflation

Cat No.	Application	Balloon OD	Balloon Length	Total length	Channel ID
EBD-M106230	Billiary/Esophagus/ Pylorus/Colon	6-7-8 mm	55 mm	230 cm	>2.8 mm
EBD-M108230	Billiary/Esophagus/ Pylorus/Colon	8-9-10 mm	55 mm	230 cm	>2.8 mm
EBD-M110230	Esophagus/Pylorus/ Colon	10-11-12 mm	80 mm	230 cm	>3.2 mm
EBD-M112230	Esophagus/Pylorus/ Colon	12-13.5-15 mm	80 mm	230 cm	>3.2 mm
EBD-M115230	Esophagus/Pylorus/ Colon	15-16.5-18 mm	80 mm	230 cm	>3.2 mm
EBD-M118230	Esophagus/Pylorus/ Colon	18-19-20 mm	80 mm	230 cm	>3.2 mm
EBD-0430080	Achalasia	30mm	80mm	90cm	N/A
EBD-0435080	Achalasia	35mm	80mm	90cm	N/A
EBD-0440080	Achalasia	40mm	80mm	90cm	N/A

*Non- Multi Stage Balloons


EUS & EBUS FNA SYSTEM

PRODUCT FEATURES

Superior Needle Visibility

Sharp Needle Atraumatic
Puncture

Large Needle opening pro?le
ensure better specimen



Cat No.	Needle G	Needle length	OD	Length	For Scope
EUS-0119142	19 G	0-80 mm adjustable	1.8 mm	142cm	Olympus/Pentax/Fuji EUS
EUS-0122142	22G	0-80 mm adjustable	1.6mm	142cm	Olympus/Pentax/Fuji EUS
EUS-0222073	22G	0-50 mm adjustable	1.6mm	73 cm	Olympus EBUS
EUS-0322073	22G	0-50 mm adjustable	1.6mm	73cm	Pentax EBUS

BIOPSY GUN

PRODUCT FEATURES

Trucut Type Biopsy Sample

Automatic Biopsy Sampling

Marking Every 1 CM

Easy Gripe Handle

AUTOMATIC BIOPSY SYSTEM



Cat No.	Biopsy Size	OD	Length
EBN - 0114016	20 mm x 0.8 mm	14 G	16 cm
EBN - 0114020	20 mm x 0.8 mm	14 G	20 cm
EBN - 0116016	20 mm x 0.8 mm	16 G	16cm
EBN - 0116020	20 mm x 0.8 mm	16 G	20 cm
EBN - 0118016	20 mm x 0.8 mm	18 G	16 cm
EBN - 0118020	20 mm x 0.8 mm	18 G	20 cm

ESOPHAGEAL DILATOR SET

PRODUCT FEATURES

Set of Eight & Sixteen Dilators

Total Length 85 cm

Marking Every 1 cm

Built - in Radiopaque Metal Ring

Highly Smooth and Flexible



Cat No.	Set Includes	Dilator Sizes
EDS - 0108085	Eight Dilators + One Spring Tip Wire	5,7,9,11,13,15 , 17,19 mm x 85 cm
EDS - 0116085	Sixteen Dilators + Spring Tip Wire	5,6,7,8,9,10,11,12,13,14,15,16, 17,18 ,19,20 mm x 85 cm

PALLIATOR

ESOPHAGEAL STENT

PRODUCT FEATURES

Nitinol Stent Highly Radiopaque

Covered From Inside

Easy Retrieval Using Drawstring


Dynamic Radial Force for Maximum Lumen Size

Dumbbell Shape for Anti-Migration

Flexible Construction to Conform to Anatomy

X-Ray Visible Marks

Scaled Marking on Delivery System



Cat No.	Stent OD	Flare OD	Stent Length	Delivery size
EMS-0118008	18cm	22 mm	8 cm	7 mm x 75 cm
EMS-0118010	18cm	22 mm	10 cm	7 mm x 75 cm
EMS-0118012	18cm	22 mm	12 cm	7 mm x 75 cm
EMS-0118015	18cm	22 mm	15cm	7 mm x 75 cm
EMS-0118010(AR)	18cm	22 mm	10 cm	7 mm x 75 cm
EMS-0118012(AR)	18cm	22 mm	12 cm	7 mm x 75 cm
EMS-0118015(AR)	18cm	22 mm	15 cm	7 mm x 75 cm

*Other sizes are available on request

PALLIATOR

BILIARY STENTS

PRODUCT FEATURES

Nitinol Stent

Small cell size


Dynamic Radial Force

Braided Tubing of Delivery System

Flexible Construction to Conform Anatomy

Zero Shortening at Proximal end

Three Marking on Delivery System for
Reliable Positioning





Cat No.	Stent Type	stent OD	Stent Length	Delivery size
EMS-0210006	Uncovered	10 mm	6 cm	8Fr x 180cm
EMS-0210007	Uncovered	10 mm	7cm	8Fr x 180cm
EMS-0210008	Uncovered	10 mm	8 cm	8Fr x 180cm
EMS-0210009	Uncovered	10 mm	9 cm	8Fr x 180cm
EMS-0210010	Uncovered	10 mm	10 cm	8Fr x 180cm
EMS-0210011	Uncovered	10 mm	11 cm	8Fr x 180cm
EMS-0210012	Uncovered	10 mm	12 cm	8Fr x 180cm
EMS-0310006	Covered	10 mm	6 cm	10Fr x 180cm
EMS-0310007	Covered	10 mm	7 cm	10Fr x 180cm
EMS-0310008	Covered	10 mm	8 cm	10Fr x 180cm

PALLIATOR

Pseudocyst Stent

PRODUCT FEATURES

- Nitinol Stent
- Fully Covered
- With Lasso
- Braided Tubing of Delivery System
- Flexible Construction to Conform Anatomy
- Three Marking on Delivery System for Reliable Positioning



Cat No.	Stent Type	stent OD	Stent Length	Delivery size
EMS-0414020	Covered	14 mm	2 cm	10Fr x 180cm
EMS-0414030	Covered	14 mm	3 cm	10Fr x 180cm
EMS-0414040	Covered	14 mm	4 cm	10Fr x 180cm
EMS-0416020	Covered	16 mm	2 cm	10Fr x 180cm
EMS-0416030	Covered	16 mm	3 cm	10Fr x 180cm
EMS-0416040	Covered	16 mm	4 cm	10Fr x 180cm

PALLIATOR

DUODENAL / PYLORIC / COLONIC STENT

PRODUCT FEATURES

Nitinol Stent

Small Cell Size


Zero Shortening at Proximal End

Braided Tubing of Delivery System

Flexible Construction to Conform Anatomy

With Lasso

Three Marking on Delivery System for
Reliable Positioning



Cat No.	Stent Type	stent OD	Stent Length	Delivery size
EMS-0522008	Uncovered	22 mm	8 cm	10Fr x 220cm
EMS-0522010	Uncovered	22 mm	10cm	10Fr x 220cm
EMS-0522012	Uncovered	22 mm	12 cm	10Fr x 220cm
EMS-0524008	Uncovered	24 mm	8 cm	10Fr x 220cm
EMS-0524010	Uncovered	24 mm	10 cm	10Fr x 220cm
EMS-0524012	Uncovered	24 mm	12 cm	10Fr x 220cm
EMS -0622008	Covered	22 mm	8 cm	10Fr x 220cm
EMS -0622010	Covered	22 mm	10 cm	10Fr x 220cm
EMS-0622012	Covered	22 mm	12 cm	10Fr x 220cm
EMS-0310008	Covered	24 mm	8 cm	10Fr x 220cm
EMS-0310010	Covered	24 mm	10 cm	10Fr x 220cm
EMS-0310012	Covered	24 mm	12 cm	10Fr x 220cm

TRACK-WIRE GUIDE-WIRE

PRODUCT FEATURES


Nitinol Core

Radiopaque Distal-Tungsten Filled Tip (5 cm)

PTFE Jacket

Hydrophilic Coated

Twin Color



Cat No.	Diameter	Length
EGW- 0225450	0.025	450 cm
EGW- 0235450	0.035	450 cm

*Other sizes are available on request

STARK-WIRE GUIDE-WIRE

PRODUCT FEATURES


Nitinol Core

Highly Radiopaque Distal 5 cm Tip

PTFE Jacket

Hydrophilic Coated

Twin Color & Highly Lubricious



Cat No.	Diameter	Length
EGW- 0125450	0.025	450 cm
EGW- 0135450	0.035	450 cm

*Other sizes are available on request

SPHINCTEROTOMES

PRODUCT FEATURES

Monofilament Cutting Wire


Triple Lumen Version

Color Marking

Separate Port for Contrast & Guidewire

Pre-curved

Ergonomical Shape Handle



Cat No.	OD	Cutting wire	Type
ESP- 0124020	2.4 mm	20 mm	Standard
ESP- 0124030	2.4 mm	30 mm	Standard
ESP- 0224007	2.4 mm	7 mm	Needle Knife
ESP- 0324020	2.4 mm	20 mm	Insulated Wire
ESP- 0324030	2.4 mm	30 mm	Insulated Wire

BILIARY DILATOR

PRODUCT FEATURES

Tapered and Rounded Ends

Well Balanced Wall Thickness

Fully Radiopaque



Cat No.	Tip OD	Dilator OD	Dilator Length
EBD- 0107200	4 Fr	7 Fr	200 cm
EBD- 0185200	5.5 Fr	8.5 Fr	200 cm
EBD- 0110200	6 Fr	10 Fr	200 cm

DRAIN - X

DRAINAGE CATHETER - INTRODUCERS - NBD CATHETER

PRODUCT FEATURES

Various Type and Sizes Available



Tapered and Rounded Ends

Well Balanced Wall Thickness

Radiopaque Material

Flaps or Pigtails

Single Step 8.5 Fr & 10 Fr Introducer



Catalog No.						Size
Biliary Amsterdam/ Curved Type	Biliary Double Pigtail	Pancreatic Curved	Biliary Single Pigtail	Biliary Multiflap/ Tannenbaum	Pancreatic Single Pigtail	
		EBS - 0305005			EBS - 0605005	5 Fr X 5 cm
		EBS - 0305007			EBS - 0605007	5 Fr X 7 cm
		EBS - 0305010			EBS - 0605010	5 Fr X 10 cm
EBS - 0107005	EBS - 0207005	EBS - 0307005	EBS - 0407005		EBS - 0607005	7 Fr X 5 cm
EBS - 0107007	EBS - 0207007	EBS - 0307007	EBS - 0407007		EBS - 0607007	7 Fr X 7 cm
EBS - 0107010	EBS - 0207010	EBS - 0307010	EBS - 0407010		EBS - 0607010	7 Fr X 10 cm
EBS - 0107012	EBS - 0207012		EBS - 0407012			7 Fr X 12 cm
EBS - 0107015	EBS - 0207015		EBS - 0407015			7 Fr X 15 cm
EBS - 0185005	EBS - 0285005		EBS - 0485005	EBS - 0585005		8.5 Fr X 5 cm
EBS - 0185007	EBS - 0285007		EBS - 0485007	EBS - 0585007		8.5 Fr X 7cm
EBS - 0185010	EBS - 0285010		EBS - 0485010	EBS - 0585010		8.5 Fr X 10cm
EBS - 0185012	EBS - 0285012		EBS - 0485012	EBS - 0585012		8.5 Fr X 12cm
EBS - 0185015	EBS - 0285015		EBS - 0485015	EBS - 0585015		8.5 Fr X 15cm
EBS - 0110005	EBS - 0210005		EBS - 0410005	EBS - 0510005		10 Fr X 5cm
EBS - 0110007	EBS - 0210007		EBS - 0410007	EBS - 0510007		10 Fr X 7 cm
EBS - 0110010	EBS - 0210010		EBS - 0410010	EBS - 0510010		10 Fr X 10 cm
EBS - 0110012	EBS - 0210012		EBS - 0410012	EBS - 0510012		10 Fr X 12 cm
EBS - 0110015	EBS - 0210015		EBS - 0410015	EBS - 0510015		10 Fr X 15 cm
EBS - 0107250	NBD Catheter Pigtail with loop					7 Fr X 250 cm
EBS - 1005220	Pusher					5 Fr X 220 cm
EBS - 1007220	Pusher					7 Fr X 220 cm
EBS - 1008220	Introducer Set					8.5 Fr X 220 cm
EBS - 1010220	Introducer Set					10 Fr X 220 cm

*Other length e.g. 4 cm, 6cm, 8cm, 9cm, 11cm, 13cm can be made on order

**All sizes are available in sets also. Please add “S” after Catalogue Nos “e.g. : EBS - XXXXXXXX-S


STONE EXTRACTION BASKET

PRODUCT FEATURES

Four Wire Construction

Braided Wires of Hard Steel

separate Flushing / Contrast port



Cat No.	Sheath OD	Basket OD	no. of Wires
EBS- 0224025	2.4 mm	25 mm	Four wire
EBS- 0224030	2.4 mm	30 mm	Four wire
EBS- 0218020	1.8 mm	20 mm	Four wire

MSE-MULTI STAGE EXTRACTOR STONE EXTRACTION BALLOON

PRODUCT FEATURES

Triple Lumen Variant

Multi Diameter Balloon Inflation

X-Ray Visible Marks

Very Strong Balloon Material



Cat No.	Lumen	Sheath OD	Injection	Inflated Diameter
EEB- 0123100	Triple Lumen	2.8 mm	Above	11.5mm - 16mm
EEB- 0223200	Triple Lumen	2.8 mm	Below	11.5mm - 16mm

LITHOTRIPTER - MANUAL


PRODUCT FEATURES

Light Weight

Easy Rotation

Compatable with all Types of Basket

Good for Impacted Stone



Cat No.	Discription	Sheath OD	Length
EML-0160090	Handle with Spiral	6 mm	90 cm
EML-0325200	4 Wire Basket	2.3 mm	200 cm
EML-0330200	4 Wire Basket	2.3 mm	200 cm
EML-0560090	Metal Spiral	6 mm	90 cm

MSE-MULTI STAGE EXTRACTOR STONE EXTRACTION BALLOON

PRODUCT FEATURES

Triple Lumen Variant

Multi Diameter Balloon Inflation

X-Ray Visible Marks

Very Strong Balloon Material



Cat No.	Discription	OD	Length
EML-0237210	Handle with Spiral	3.7 mm	210 cm
EML-0425210	4 Wire Basket	25 mm	210 cm
EML-0430210	4 Wire Basket	30 mm	210 cm
EML-0537210	Metal Spiral	3.7 mm	210 cm
EML-0625210	Metal Spiral with Basket	3.7 mm/25mm	210 cm
EML-0630210	Metal Spiral with Basket	3.7 mm/30mm	210 cm

P.E.G

PERCUTANEOUS ENDOSCOPIC GASTROSTOMY FEEDING TUBE

PRODUCT FEATURES

Pull Technique

Thin Wall Tubing

Large Inner Lumen



Cat No.	OD	Kit Contents
EPG-0120030	20 fr	Mini Set -PEG Feeding tube, Feeding Adapter, Retention Ring, Scalpel, Introducer, Tubing Clamp, Looped Placement Wire
EPG-0124030	24 fr	Mini Set -PEG Feeding tube, Feeding Adapter, Retention Ring, Scalpel, Introducer, Tubing Clamp, Looped Placement Wire
EPG-0220030	20 fr	Full Set - PEG Feeding tube, Feeding Adapter, Retention Ring, Scalpel, Introducer, Tubing Clamp, Looped Placement Wire, Snare, Scissors, Hemostat, Guaze, Drape, Needle, Syringe, Povidon Iodine,
EPG-0224030	24fr	Full Set - PEG Feeding tube, Feeding Adapter, Retention Ring, Scalpel, Introducer, Tubing Clamp, Looped Placement Wire, Snare, Scissors, Hemostat, Guaze, Drape, Needle, Syringe, Povidon Iodine,

Gastrostomy Feeding Tube

PRODUCT FEATURES

Medical Grade %100 Silicon Tube

Thin Wall Tubing

Large Inner Lumen

Non-traumatic Distal Tip



Cat No.	Tube OD	Tube Length	Balloon Volume
EPG-2112024	12 Fr	24 cm	5 ml
EPG-2114024	14 Fr	24 cm	5 ml
EPG-2116024	16 Fr	24 cm	20 ml
EPG-2118024	18 Fr	24 cm	20 ml
EPG-2120024	20 Fr	24 cm	20 ml
EPG-2122024	22 Fr	24 cm	20 ml
EPG-2124024	24 Fr	24 cm	20 ml



Pneumatic Inflation Device



Channel Cap with Flushing Tube




Disposable Cytology Brush



Scope Cleaning Brush



Endoscope Bite Guard



Balloon Inflation Device



Viza med



sales@viza-med.com



+963 956 079 779

