

Laparoscopy accessories

Trocar & Reducer&Cap



SK502-401

#code :

Specification:

Φ5.5X 100mm/Φ10.5mmX100mm Φ12.5X100mm, Φ15X 100mm
(Flip type)




SK401-560/561/562

#code :

Specification:

Φ10.5-Φ5.5mm, 120mm/100mm/22mm




SK401-504

#code :

Specification:


Φ5.5X 100mm/Φ10.5mmX100mm Φ12.5X100mm, Φ15X 100mm
(Magnetic type)

Trocar & Reducer&Cap




#code : SK401-508/555

Specification: $\Phi 10.5-\Phi 5.5\text{mm}$



#code : SK401-506


Specification: $\Phi 5.5 \times 100\text{mm} / \Phi 10.5\text{mm} \times 100\text{mm} \ \Phi 12.5 \times 100\text{mm}, \ \Phi 15 \times 100\text{mm}$
(Titanium alloy type)



#code : SK401-551/551A

Specification: $\Phi 10.5 \& \ \Phi 5\text{mm},$

Monopolar HF Electrode



#code : SK401-130/131/135/130A/132/136/133/134/137

Specification: $\Phi 5 \times 330\text{mm}$



#code : SK401-138/138A/139/139A

Specification: Suction irrigation $\Phi 5 \times 330\text{mm}$



#code : SK401-138/138A/139/139A

Specification: Cable for Monopolar $\Phi 4 \times 3000\text{mm}$



#code : SK401-577/577A

Specification: ordinary connector 6mm/8mm

Suction& Irrigation, Clip Applier & Bipolar



#code : SK401-140

Specification: Suction&Irrigation Φ5 /Φ10 X 330mm
Push type



#code : SK401-141

Specification: Suction&Irrigation Φ5 /Φ10 X 330mm
Spring type



#code : SK401-107

Specification: Polymer Clip Applier Φ5 /Φ10 X 330mm



#code : SK401-105

Specification: Polymer Clip Applier Φ5 /Φ10 X 330mm with locking

Trocar & Reducer&Cap



#code :

SK401-115/115A/117

Specification:

Ø5 X 330mm, single action/double action, Dissect/Fenestrated



#code :

SK401-119

Specification:

Ø4 X 3000mm, Bipolar cable



#code :

SK401-101/102/103

Specification:

Titanium Clip Applier

Needle-holding forceps & Hem-o-Lok clips



V 型持针钳 V shape Needle Holder Φ 5 X 330mm
SK401-122 直头 Straight
SK401-122A 左弯 Curved left
SK401-122B 右弯 Curved right
SK401-122C 自动复位 Self-righting
SK401-122D 复位夹线 Self-righting with teeth

#code : SK401-122/A/B/C/D

Specification: Needle-holding forceps Φ 5 X 330mm Material:Stainless Steel



V 型持针钳 V shape Needle Holder Φ 5 X 330mm
SK401-123 直头 Straight
SK401-123A 左弯 Curved left
SK401-123B 右弯 Curved right
SK401-123C 自动复位 Self-righting
SK401-123D 复位夹线 Self-righting with teeth

#code : SK401-123/A/B/C/D

Specification: Needle-holding forceps Φ 5 X 330mm Material:Stainless Steel



SK401-127 A型持针钳(钛合金) A shape(Titanium Alloy) Needle Holder Φ 5 X 330mm
SK401-127 直头 Straight
SK401-127A 左弯 Curved left
SK401-127B 右弯 Curved right
SK401-127C 自动复位 Self-righting
SK401-127D 复位夹线 Self-righting with teeth

#code : SK401-127

Specification: Needle-holding forceps Φ 5 X 330mm Material:Titanium Alloy



SK401-127 A型持针钳(钛合金) A shape(Titanium Alloy) Needle Holder Φ 5 X 330mm
SK401-127 直头 Straight
SK401-127A 左弯 Curved left
SK401-127B 右弯 Curved right
SK401-127C 自动复位 Self-righting
SK401-127D 复位夹线 Self-righting with teeth

#code : SK401-127A/B/C/D

Specification: Needle-holding forceps Φ 5 X 330mm Material:Titanium Alloy

Needle-holding forceps & Hem-o-Lok clips




#code : SK401-125/A/B/C/D

Specification: Φ 5.5 X 330mm Material:Stainless Steel



#code : SK401-126/A/B/C/D

Specification: Φ 5.5 X 330mm Material:Stainless Steel



#code : SK401-128

Specification: Φ 5.5 X 330mm Material:Titanium Alloy



#code : SK401-128A/B/C/D

Specification: Φ 5.5 X 330mm Material:Titanium Alloy

Needle-holding forceps & Hem-o-Lok clips



#code : GCJZ-XL*6

size: XL

Specification: 6 clips in one cartridge,
20 cartridges in one box.



#code : GCJZ-L*6

size: L

Specification: 6 clips in one cartridge,
20 cartridges in one box.



#code : GCJZ-ML*6

size: ML

Specification: 6 clips in one cartridge,
20 cartridges in one box.

Knot Pusher & Fan-shaped forceps



#code : SK401-149

Specification: $\Phi 5 \times 330\text{mm}$
(steel head)



#code : SK401-150

Specification: $\Phi 5 \times 330\text{mm}$
(plastic head)



#code : SK401-152

Specification: $\Phi 5 \times 330\text{mm}$
(slot type)

Knot Pusher & Fan-shaped forceps



#code : SK401-152

Specification: $\Phi 5 \times 330\text{mm}$
(slot type)



#code : SK401-212

Specification: $\Phi 5 \times 330\text{mm}$
(three leaves)



#code : SK401-213

Specification: $\Phi 5 \times 330\text{mm}$
(three leaves)

mm Laparoscopic Forceps, Graspers 5mm



#code : SK401-161

Specification: Ø5mm scissors, straight with single action



#code : SK401-184

Specification: Ø5mm Alligator grasping forceps 14mm



#code : SK401-166

Specification: Ø5mm scissors, curved with double action



#code : SK401-186

Specification: Ø5mm Maxi grasper

mm Laparoscopic Forceps, Graspers 5mm



#code : SK401-167

Specification: Ø5mm scissor hook



#code : SK401-187

Specification: Ø5mm Babcock forceps



#code : SK401-172

Specification: Ø5mm straight dissecting forceps



#code : SK401-189

Specification: Ø5mm Oviduct grasping forceps

5mm Laparoscopic Forceps, Graspers 5mm



#code : SK401-163

Specification: Φ 5mm scissors, straight with double action



#code : SK401-191

Specification: Clinch grasping forceps



#code : SK401-162

Specification: Φ 5mm scissors,
single-action with ball



#code : SK401-192

Specification: Φ 5mm
fenestrated grasping forceps

mm Laparoscopic Forceps, Graspers 5mm



#code : SK401-170

Specification: Φ 5mm Biopsy forceps, spoon jaw



#code : SK401-193

Specification: Φ 5mm
Blunt grasping forcep



#code : SK401-173A

Specification: Φ 5mm Maryland 12mm



#code : SK401-195A

Specification: Φ 5mm
Fenestrated forceps 20mm

mm Laparoscopic Forceps, Graspers 5mm



#code : SK401-173

Specification: Φ 5mm Maryland 17mm



#code : SK401-195B

Specification: Φ 5mm
SK401-195B



#code : SK401-173B

Specification: Φ 5mm Maryland 14mm



#code : SK401-200

Specification: Φ 5mm
knot forceps single action

mm Laparoscopic Forceps, Graspers 5mm



#code : SK401-173C

Specification: Φ 5mm Maryland 18mm



#code : SK401-201

Specification: Φ 5mm
B-knot forceps



#code : SK401-173D

Specification: Φ 5mm Maryland 20mm



#code : SK401-202

Specification: Φ 5mm
S-knot forceps

mm Laparoscopic Forceps, Graspers 5mm



#code : SK401-173E

Specification: Φ 5mm Maryland forceps,
25mm



#code : SK401-173F

Specification: Φ 5mm Maryland forceps
16mm blunt



#code : SK401-173F

Specification: Φ 5mm Maryland forceps
14mm



#code : SK401-173G

Specification: Φ 5mm Maryland forceps
14mm

mm Laparoscopic Forceps, Graspers 5mm



#code : SK401-179

Specification: Φ 5mm Cobra grasping forceps double action with 2+3 teeth



#code : SK401-180

Specification: Φ 5mm Cobra grasping forceps double action with 2+3 teeth



#code : SK401-181

Specification: Φ 5mm Grasping forceps 19mm single action with 2+3



#code : SK401-183

Specification: Φ 5mm Grasping forceps 12mm single action with 1+2

10mm Laparoscopic Forceps, Graspers 10mm



#code : SK401-304

Specification: Φ 10mm
large straight scissor, double action



#code : SK401-305

Specification: Φ 10mm
Dissecting forceps



#code : SK401-306

Specification: Φ 10mm
(grasping forceps)



#code : SK401-310

Specification: Φ 10mm
(Lithotomy forceps)

10mm Laparoscopic Forceps, Graspers 10mm



#code : SK401-302

Specification: $\Phi 10\text{mm}$ large straight scissor




#code : SK401-303

Specification: $\Phi 10\text{mm}$ large straight scissor, single action



#code : 804.043

Specification: $\Phi 2.8 \times 280\text{mm}$
(punch grasping forceps)




#code : SK504


Specification: Uterus - cutter and
allessories one set

Disposable Electric Endoscopic Linear Cutting Staplers and reload

Code : EECRA-45B
EECRA-45G
EECRA-45H
EECRA-45W
EECRA-60B
EECRA-60G
EECRA-60H
EECRA-60W



Code : EECRB-45B
EECRB-45G
EECRB-45H
EECRB-45W
EECRB-60B
EECRB-60G
EECRB-60H
EECRB-60W



Code : EECSA250-45Y
EECSA250-60Y
EECSA250-45P
EECSA250-60P
EECSA350-45P
EECSA350-60P



Disposable Trocar for Laparoscope (Bladeless)

Code : GCCQA-5
 GCCQA-10
 GCCQA-12
 GCCQA-15
 GCCQB-5
 GCCQB-10
 GCCQB-12



Disposable Non-absorbable Closure Clip (Polymer)

Code : GCJZ-ML*2
 GCJZ-ML*4
 GCJZ-ML*6
 GCJZ-L*2
 GCJZ-L*4
 GCJZ-L*6
 GCJZ-XL*2
 GCJZ-XL*4
 GCJZ-XL*6



Disposable Circular Staplers

Code : GWA-23H
 GWA-26H
 GWA-29H
 GWA-32H



Disposable Circular Staplers for Hemorrhoids

Code : GCGZA-32
 GCGZA-34



Disposable Endoscopic Cutting Stapler and Reload

Code : G CJQB-160
G CJQB-260



Code : G CJQW-45H
G CJQW-45L
G CJQW-45M
G CJQW-60L
G CJQW-60M



Disposable Linear Cutting Staplers and Reload

Code : G CLCA-75M



Code : G CLCA-ZJ75M



Disposable Skin stapler and staple remover

Code : G CPFB-35W



Code : G CPC-D



Disposable Wound Protector

Code : QKA 60/70-60/150
 QKA 80/90-80/150
 QKA 120/130-120/250
 QKA 150/160-150/250
 QKA 180/190-180/250
 QKA 220/230-220/250



Code : QKB 40/30-15/25
 QKB 50/40-25/25
 QKB 70/70-35/25
 QKB 70/70-35/35



Disposable Retrieval Bag

Code : QWA-100
 QWA-130
 QWA-80



Disposable Anal Hemorrhoid Ligation Device (with Elastic thread)

Code : TZ-3/4



Disposable Anal Hemorrhoid Ligation Device (with rubber band)

Code : TZ-3/4



Disposable Minimally Invasive Fascia Closure Device

Code : GCJM



Disposable Thoracentesis apparatus

Code : GCXCA-6/10/12




Product Description: *Manual Control Forceps Vessel Sealer*

Specification & Model: **1201**

Jaw Size mm: **≥18**

Coagulation length mm: **16.5**

Cutting Length mm: **14**



Product Description: *Manual Control 13mm Vessel Sealer*

Specification & Model: **4218**

Shaft Length mm: **180**

Shaft Diameter mm: **13**

Jaw Size mm: **≥18**

Coagulation length mm: **36**

Cutting Length mm: **34**



REF: **TCLT300**

Size: **Large**

Colour: **Green**

PCS/Cartridge: **6**

Cartridge / Box: **20**





Product Description: **Manual Control 5mm Vessel Sealer**

Specification & Model: **LB05**


Shaft Length mm: **370/440**

Shaft Diameter mm: **5**

Jaw Size mm: **≥18**

Coagulation length mm: **20.5**

Cutting Length mm: **18.5**



Product Description: **Manual Control 5mm Vessel Sealer**

Specification & Model: **LXB05**

Shaft Length mm: **370/440**

Shaft Diameter mm: **5**

Jaw Size mm: **≥18**

Coagulation length mm: **20.5**

Cutting Length mm: **18.5**





viza-med
VISION & PRECISION